

Peak identification and alignment

Yutaka Yasui, Ph.D.

Dale McLerran, M.S.

**Division of Public Health Sciences
Fred Hutchinson Cancer Research Center**

Collaborators

NCI Early Detection Research Network

Sudhir Srivastava

Eastern Virginia Medical School

John Semmes

George Wright, Jr.

Bao-Ling Adam

Fred Hutchinson Cancer Research Center

Bree Mitchell

Phil Gafken

Paul Lampe

Johanna Lampe

Dale McLerran

Margaret Pepe

Tim Randolph

Yinsheng Qu

Li Hsu

Mary Lou Thompson

Mark Thornquist


Yan Liu

Toana Kawashima

Marcy Winget

John Potter

Ziding Feng



Recognition of a problem
in the current approach





Solution of the problem
via
modifications / new developments



Mass-spectrometry technology (MALDI, SELDI)


Matrix (Surface Enhancement = SELDI)




Mass/Charge



An example of SELDI output



Analysis Steps

- Calibration from TOF to mass/charge (M/Z)
 - Baseline subtraction / Normalization
 - Peak identification
 - Peak Alignment
 - Search for signature profiles
- 

1. Calibration

Conversion of TOF to mass/charge

$$\frac{\frac{m}{z} + pm}{\text{Voltage}} = \alpha(\text{TOF} - \beta)^2 + \gamma$$

- Measure **TOF** of 7 (or 5) peptides with known **m/z** values
- Fit the above equation and estimate the parameters (α , β , γ)
- Apply the derived equation to convert **TOF** to **m/z**



Calibration Issues

Goodness of fit with the 7 (or 5) standard peptides

Check the goodness of fit by eliminating one standard peptide: identify any “bad” standard peptide(s)

**Over the course of an experiment,
what is the optimal schedule of calibration?
(once, multiple times, everyday, ...)**

After Calibration




2. Baseline subtraction & Normalization

Subtract the amount of intensity
inflated by matrix


Scale the intensity to normalize spectra
(total ion current)

Before baseline subtraction




**Baseline intensity
due to matrix**


After baseline subtraction



After Baseline Subtraction & Normalization



~ 50,000
(m/z , BLsubtracted-
normalized intensity
points



3. Peak Identification

A mathematical definition
of peak locations


A critical issue in SELDI/MALDI-TOF analyses



Two imprecision problems

1. Imprecise measurements of **mass/charge** values (X-axis)
2. Imprecise measurements of **intensity** values (Y-axis)

Properties of SELDI / MALDI-TOF output

Shift ± 0.1 -0.2% of m/z (QC, not uniform)







1. Define peaks
(similarly to how mass spectra are read/utilized)

2. Fix miss-aligned peaks

Directed approach (case-control blinded)




Peak identification

Ask at each point:


*“Is it the highest point
in the neighborhood
of $\pm N$ points?”*

YES = a peak


NO = not a peak



Peak identification results



Peak refinement



1. Moving median smoother with a wide span to obtain typical intensity values $\{(m/z, s)\}$ locally
2. Moving median smoother with a wide span to obtain typical absolute deviation $\{(m/z, d)\}$ locally
3. Define a peak if $\text{intensity} > s + K * d$

Peak identification

- Identify peaks in each spectrum
- The number of peaks per sample is $\sim 1,000$ ($\sim 2\%$ of the original 50,000 points)
- High- and **low-intensity** peaks
(vs. ~ 150 peaks by Coombes)
- Now, align peaks across samples
(**Alignment**)

4. Peak Alignment

Correction of miss-aligned peaks
across spectra

sample A

sample B

sample C

sample D



Save as the 1st
aligned
mass/charge &
its intensity in
the aligned
dataset




$\leftrightarrow \pm P\%$ of the mid mass/charge
 \downarrow
1st aligned mass/charge, X_1



Save as the 2nd
aligned
mass/charge &
its intensity in
the aligned
dataset

2nd aligned mass/charge, X_2





**Save as the 3rd
aligned
mass/charge &
its intensity in
the aligned
dataset**

3rd aligned mass/charge, X_3



1st

2nd

3rd

...

...

The aligned dataset
for
searching
signature markers
profiles

Completion of pre-
analysis processing


Yasui et al. J. Biomed. & Biotech
(Special Issue on Proteomics) 20



A new modification of our alignment algorithm by Dale McLerran


Instead of using the number of peaks,

use a weighted sum (ws) of peak counts in a window




- If $b < \min(a, c)$ then $ws = 0$
- If $b > \max(a, c)$ then $ws = a + \lambda b + c$
where $\lambda = b / \max(a, c)$
- Otherwise $ws = a + b + c$

Original alignment method peak locations



Overweight centers




Imprecise measurements of **intensity** values (Y-axis)

Identify and quantitate sources of variations

Replicate the measurements at high-variation sources

Variance components assessment

Repeated measurements of a single QC sample



Sources of
intensity variation

V_{day} = Day-to-Day

V_{chip} = Chip-to-Chip

V_{well} = Well-to-Well

Reduction of CV by averaging 3 replicates

$$CV = (\text{Variance of intensity})^{1/2} / \text{Mean}$$

$$= (V_{\text{day}} + V_{\text{chip}} + V_{\text{well}})^{1/2} / \text{Mean}$$

If spotted on 3 wells of a chip,

$$CV_3 = (V_{\text{day}} + V_{\text{chip}} + 1/3 \times V_{\text{well}})^{1/2} / \text{Mean}$$

Reduction of CV by averaging 3 replicates

If measured on 3 different days,

$$\begin{aligned} CV_3 &= \{1/3 \times (V_{\text{day}} + V_{\text{chip}} + V_{\text{well}})\}^{1/2} / \text{Mean} \\ &= CV / \sqrt{3} \end{aligned}$$

(e.g., $CV = 20\%$ then $CV_3 = 11.5\%$)

Summary

- Proper calibration/normalization is critical
- Imprecise measurements of m/z values necessitate the identification and alignment of peaks
- Simple algorithms have been developed and available
- Further refinements and alternative approaches are possible (need quick developments even if not optimal)
- Replication alleviates the imprecision problem of intensities

